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March 11, 2012

# 1 Introduction

The **econum** R package is a collection of R functions used at the EcoNum lab to perform various tasks and calculations. It contains routines for (note that much of these function still need to be added!):

- ▷ Interface many scientific devices like ph-meters, titrators, scales, CO<sub>2</sub>-recorder, IKS computers, etc.,
- ▷ Controlling the proprietary chemostat, including display and diagnostic graphs,
- ▷ Pilot the automatic titrators for the alkalinity and calcium titrations,
- ▷ Calculate concentrations and stats on measurements done using a spectrophotometer,
- ▷ Perform various calculations and conversions related to ecophysiological studies of aquatic invertebrates,
- ▷ And much more...

## 1.1 Installation of the econum R package

Get R (at least 2.10.0) and install it. Then, get either the source (for Linux and Mac OS X) of the **econum** package ('econum\_x.y-z.tar.gz'), or the Windows binaries ('econum\_x.y-z.zip') from 'EN-Dolphin'. Install them inside R (see R online help and manuals to learn how to do that). You need administrative rights on the machine to do this install.

Under Windows (95 -> XP), if you intend to pilot the parallel interface card, you must also copy 'inpout32.dll' and 'lpttcl.dll' from '/bin' subdirectory of the **econum** package to 'c:\Windows\system32' (or whatever directory where system dlls are installed).

## 2 General instructions for working in the lab

This chapter covers general instructions that are supposed to be known for working in the lab, but that are always good to check. It covers general actions to do good laboratory work. TODO: continue this and cover good practices in the laboratory topic.

- ▷ General behavior in the lab and security issues: no food, no drink, correct clothes, always keep working area clean and in order, work under the “hotte”, know basic security issues (where is the extintor? what to do in case of ..., etc.)
- ▷ The water: various water qualities and the importance of using the right water quality.
- ▷ Cleaning equipment and glassware.
- ▷ Logging all activities in a logbook.
- ▷ What else?

### 2.1 Marine chemistry: basic knowledges on concentration in seawater

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### 3.1.3 Références

[refs TODO!].

## 3.2 Measuring zooxanthellae with the FlowCAM

TODO: make the protocole here. The following part is just a transcript of instructions received from Fluid Imaging...

### 3.2.1 Operation of the FlowCAM using a 50 µm Flow cell

For processing samples within 15-20 min using a 20x objective and 50 µm flowcell:

1. Shorten the length of the tubing between the funnel and the flow cell as much as possible (~1cm).
2. If possible use a smaller funnel with a narrow base (less overall volume – one possible idea would be to cut the top end of a disposable plastic pipet that has a very narrow stem diameter).
3. Add the sample and set the pump to PRIME.
4. Once the sample has reached the flow cell, process the sample using autotrigger mode for the necessary period of time (~5-10 mins).
5. Once the run has finished, remove excess sample in the funnel with a pipet (disposable), add and rinse funnel with filtered seawater. Have the pump set to prime. Rinse funnel with filtered seawater and continue to remove fluid from the funnel while the pump is priming with a pipet – continue for 30 seconds until FSW has entered the 50 µm depth flow cell.

### *3 Laboratory protocols*

6. Remove the FSW from the funnel to the neck of the funnel. Add your next sample and prime the pump once again until the sample has entered the flow cell (monitor using setup and focus), then begin to run the next sample.
7. With a 5-10 min run using autotrigger, and these rinsing techniques the processing time of each sample should be approximately 15-20 minutes.

## 4 Computer control of scientific devices

At EcoNum, we use scientific devices connected to computers as much as possible. On one hand, it allows great flexibility in automatizing even complex tasks by writing computer programs to pilot the scientific devices, and in the other hand, it is fast and reliable to collect data from instruments directly into the computer as much as possible. When you take note of your measurements, then introduce them manually in the computer, you waste your time, but there is also a possibility to make mistakes when you copy the data.

First few sections deal with basic interfaces management (ph-meters, controller cards, scales, etc.), and following sections cover more elaborate implementations like the automatic titration procedures, and the chemostat.

### 4.1 The parallel interface card

The parallel interface card is very simple to pilot. It is plugged to the computer with a parallel cable and it allows to master eight channels (with relays) corresponding to the eight bits of a byte. Thus, sending one byte to the card through the parallel cable instructs to open (bit is 1) or close (bit is 0) respective channels. The **econum** package contains functions to pilot that card through a special LPT driver under Tcl/Tk. This driver is only working under Windows 95 -> XP. The card must be connected to usual LPT1, LPT2, or LPT3 parallel ports (no exotic parallel ports supported yet). Its use is very simple through the unique `para()` function. Basically, you specify a vector of eight logicals indicating the state of each channel you want, the parallel port where the card is connected<sup>1</sup>, and you change the state of the card instantly. Additionally, you can specify the argument `change = FALSE` and in this case, the function does not change the state of the eight channels but it reads the current state of the card. Here is an example:

```
> options(para.port = 1)      # If your card is connected to LPT1
> para(c(T,F,F,T,F,F,F,T))  # Channels 1, 4, 8 ON, the others OFF
> para(change = FALSE)       # Read current state
> para(c(F,F,F,F,F,F,F,F))  # All channels OFF
> para(change = FALSE)       # Make sure current state is changed
```

---

<sup>1</sup>It is much more convenient to use the `options(para.port = X)` instruction before using `para()`, where X is the parallel port 1, 2 or 3. So that, you set up the port once and you do not need to specify it again and again each time you call `para()` (see the example).

# 5 Annexes

## 5.1 A couple of recommendations

- ▷ Be consistent with the spelling of your name in all your scientific literature. In particular, decide if you indicate initial al middle name(s) or not, consistently (J. Doe *versus* J.F. Doe). Changing this from paper to paper makes it difficult to follow your work and you will be listed under different authors in the online databases.
- ▷ Choose a scientific journal *before* you write your paper. Make sure to read which topics the journal covers and to write your paper accordingly. Read carefully author's instructions and follow them from the beginning.
- ▷ Always check copyright policies of the journals you submit paper, and later on, when you want to send reprints (possibly in PDF format) or publish them on Web sites. It is not always possible to do so. Sometimes, you can only publish the summary, but not the full paper. If you want to diffuse a work largely and easily (a white paper whose goal would be to be accessible by everybody, for instance), choose a publisher that allows you to freely publish it on the Web.
- ▷ Make sure all your scientific literature (from the most confidential internal note or report to published books) has page numbers, correct and complete references, author names, date, and possibly, version, ISBN, or copyright information if required.
- ▷ Usual rules for page numbering is to number *every* page (even blank ones) and start at one for the recto of first printed page in arabic numerals. Table of contents and index may be numbered separately, possibly in roman numerals to avoid confusion. When you decide how to number your pages, consider if your work will be delivered as a PDF. In this case, it is much better to decide to number all pages in a single sequence with arabic numerals, starting at one from the title page. That way, page numbering will be the same as the page number displayed in, say, Adobe Acrobat Reader or GostView. The L<sup>A</sup>T<sub>E</sub>X styles "article" and "report" do this by default, but not the "book" style, for instance.
- ▷ To distribute your work online, always prefer Adobe's Portable Document Format (PDF). Even for presentations, like those produced with Microsoft Powerpoint, do prefer PDF for distribution, but take care that complex transitions and animations are lost (it is always better to avoid them, when possible, in presentations)!

- ▷ When converting your document to PDF, always keep in mind the following points: (1) will your vector drawings be converted into raster images, and if so, how would they look like on screen and printed? (2) What resolution for images shall you choose? The optimal one is different for documents to be viewed on screen or to be printed. If you need both, consider something like 150 dpi, which is often a good compromise, but test the PDF to get an idea of the result, both on screen and after printing it on a good laser and a good inkjet printers. (3) The *shrink to fit* function in Adobe Acrobat Reader is nice for printing non-ISO paper size, and in particular to cope with little difference between standard A4 format and North American Letter format. (4) PDF uses built-in compression. Documents are often smaller than in their original format (especially if you embed large images in them), after choosing a reasonable resolution of course (see point 2). You don't need to zip, or compress the PDF with another mean, and the gain would be marginal. (5) Make sure your PDF converter creates a readable table of contents, and can manage hyperlinks (URLs, but also internal cross-references, links to titles in the document from the table of content, and links to corresponding pages in the index. Not all PDF converters take this into account, especially those using a virtual printer driver. OpenOffice converter does a good job, as well as Adobe Acrobat Distiller for Microsoft Office documents, and L<sup>A</sup>T<sub>E</sub>X when using the **hyperref** package. (6) Always check if the fonts you use are available and look for the possibility to embed most exotic fonts in your document (display the PDF on a machine that does not have that font installed). (7) Some old PDF converters, like Ghostscript < 6.5 lack support for Type1 fonts and transform exotic fonts into ugly raster representations. Take care of this and consider a different converter if you notice an obvious degradation of the font in the final PDF document. (8) Does your PDF converter and reader manage to get good-looking anti-aliased images, especially for your vector drawings? Not all couples of converters/readers do this. If you have a reader that does not anti-alias text and images, try to find one that does it: it dramatically increases the readability of your document, especially on low-resolution screens. (9) Look at metadata (title, author, keywords, etc.) embedded in your PDF document. Are these correct and up-to-date? If not, look in the documentation of your PDF converter to determine how you can change this (but it usually takes the data from the *properties* defined for the initial document. (10) Use appropriate format for all your pictures. A vector graphic format (SVG, WMF/EMF, PDF, ...) is better for diagrams and graphs. A raster format should be used for the others. There are two types of raster formats: those who compress data without loss of information (GIF, PNG, ...) that are most suitable for screenshots or any images containing text and large areas of uniform colors with sharp edges, and those with lossy compression (typically JPEG) that do an excellent job with maximum compression for photos that have to be viewed unmodified by humans. Even for photos, do not use JPEG if the image must be modified, or some measurements need to be done on them using an image analysis software. Check the number of colors or gray levels the format can handle. In particular, avoid GIF format if

you need more than 256 different colors in your image. Take care that many distillers use a JPEG compression for *all* images in your document. There is a little script [here](#) that deactivates JPEG compression for selected images in the process of creating a PDF.

- ▷ You can test your PDF converter/distiller and reader with the [testflow](#) document if you use T<sub>E</sub>X. That document contains a series of commands and features that are known to cause problems in the production of a PDF version of your L<sup>A</sup>T<sub>E</sub>X work.
- ▷ Always think at your production chain and consider solutions that are *traceable* and *reproducible*. Here are a couple of recommendations.